

# Wilma LabGlass Economy NMR Tube

**R**ecently some of our competitors altered our MHz rating standard, especially for High-Throughput Economy Tube, by advertising inferior product for higher frequency experiments.

To avoid lowering the SNR, we offer this complete **NMR Economy Tube Cross Reference Guide**. From this head-to-head comparison, Wilma LabGlass proves its 60 year philosophy again - we are committed to providing customers with cost-effective NMR consumables of the highest quality.



MHz Rating	Outer Diameter (mm)	Wilma LabGlass <sup>1</sup>		NORELL <sup>®</sup>		New Era		Kontes <sup>®</sup>	
		Model Number	Price <sup>2</sup>	Model Number	Price <sup>5</sup>	Model Number	Price <sup>5</sup>	Model Number	Price <sup>5</sup>
700	5.0	WG-1242-7	\$14.20						
700	5.0	WG-1242-8	\$17.12						
600	5.0	WG-1241-7	\$6.76			NE-SL5-7	\$10.50		
600	5.0	WG-1241-8	\$8.08			NE-SL5-8	\$11.93		
500	5.0	WG-1235-7	\$6.13	509-UP-7	\$16.92	NE-UL5-7	\$8.92	897150-0007	\$7.29
500	5.0	WG-1235-8	\$7.44	509-UP-8	\$19.22	NE-UL5-8	\$10.24	897150-0008	\$8.82
400	5.0	WG-1228-7	\$4.29	508-UP-7	\$8.10			897140-0007	\$5.17
400	5.0	WG-1228-8	\$5.49	508-UP-8	\$9.61			897140-0008	\$6.55
300	5.0	WG-1226-7	\$3.78	506-P-7	\$5.02	NE-HL5-7	\$5.07	897130-0007	\$4.53
300	5.0	WG-1226-8	\$5.01	506-P-8	\$6.17	NE-HL5-8	\$6.28	897130-0008	\$5.68
200	5.0	WG-1208-7	\$2.92	XR-55-7	\$4.23	NE-ML5-7	\$4.33	897120-0007	\$3.45
200	5.0	WG-1208-8	\$3.40	XR-55-8	\$5.02	NE-ML5-8	\$4.54	897120-0008	\$4.14
100	5.0	WG-1206-7	\$2.64	505-P-7	\$3.87	NE-LL5-7	\$4.01	897110-0007	\$3.15
100	5.0	WG-1206-8	\$3.09	505-P-8	\$4.59	NE-LL5-8	\$4.33	897110-0008	\$3.68
100	5.0	WG-5MM-Economy-7	\$1.72	502-7 552-7	\$1.88			897193-0000	\$3.30
100	5.0	WG-5MM-Economy-8	\$2.27	502-8 552-8	\$2.35			897193-0008	\$3.59
HT <sup>3</sup>	5.0	WG-1000-7	\$108.37	ST500-7	\$147.44	NE-RG5-7	\$121.25	897193-0050 897193-7100	\$146.49
HT <sup>3</sup>	5.0	WG-1000-8	\$134.78	ST500-8	\$174.24			897193-8050 897193-8100	\$154.14
HT <sup>3</sup>	10.0	WG-4000-7	\$333.78	1001-7	\$351.36				
TD <sup>4</sup>	10.0	WG-4001-7	\$106.99					897300-0000	\$113.06
HT <sup>3</sup>	3.0	WG-3000-7-50	\$103.99	S-3-HT-7	\$150.60				
HT <sup>3</sup>	3.0	WG-3000-8-50	\$114.59	S-3-HT-8	\$165.64				

1: Wilma LabGlass's Economy Tube is 30% more robust than any competitors' products by increasing wall thickness from 0.38 mm to 0.43 mm.

2: The price reflects 2015 US market suggested retail price. The lowest price is marked in BLUE.

3: HT refers to High-Throughput and is bulk packed.

4: TD refers to time domain NMR. TD NMR tubes are flat bottomed, and bulk packed.

5: The competitor's price is based on the same quantity as Wilma LabGlass's corresponding product.

NORELL<sup>®</sup> is a registered trademark of Norell, Inc. Kontes<sup>®</sup> is a registered trademark of Kimble Chase Life Science and Research Products LLC.



Toll Free: **1-800-220-5171**  
Intl.: **856-691-3200**



Fax

Toll Free: **1-800-220-1081**  
Intl.: **856-691-6206**



Web

[www.wilmad-labglass.com](http://www.wilmad-labglass.com)  
[cs@wilmad-labglass.com](mailto:cs@wilmad-labglass.com)