

# Lab Liquid Mixing Equipment

## • LG-MX Vortex Mixer



LG-MX Technical Data	
Power (W)	60
Shaking Movement	orbital
Motor Input (W)	58
Motor Output (W)	10
Speed Range (rpm)	0-2500
Switch	Continuous/Touch
Dimension [D×W×H] (mm)	127×130×160

- Easily Set Speed from 0 to 2500 rpm
- Select from Continuous On/Off/Touch On
- Water-proof Switch Cover
- Standard Anticorrosive Silicone Base
- Exchangeable Vacuum Feet Included for Improved Stability

Product Number	Included Items	Price (2013)
LG-MX	<ul style="list-style-type: none"> <li>✓ Vortex Mixer with Standard Silicone Base</li> <li>✓ Exchangeable Vacuum Feet (4×)</li> </ul>	<div style="border: 2px solid black; border-radius: 50%; padding: 5px; display: inline-block;"> <small>ONLY</small>  <b>\$159.14</b> </div>



Toll Free: **1-800-220-5171**  
Intl.: **856-691-3200**



Toll Free: **1-800-220-1081**  
Intl.: **856-691-6206**



[www.wilmad-labglass.com](http://www.wilmad-labglass.com)  
[cs@wilmad-labglass.com](mailto:cs@wilmad-labglass.com)