

# CELL CULTURE

## Proculture® Glass Spinner Flasks Parts & Accessories



1172 NW Blvd. Vineland, NJ 08360 USA  
Phone: 800-220-5171 / 856-691-3200 • Fax: 856-691-6206 / 800-220-1081  
[www.Wilmad-LabGlass.com](http://www.Wilmad-LabGlass.com)



*An SP Industries Company*

*Committed to Product Excellence... and Continuous Support*



***OUR CUSTOM CAPABILITIES ARE VARIED, INCLUDING CENTER NECK PORTS, SIZE & CONFIGURATION OF SIDE NECKS & CUSTOM INTERNAL COMPONENTS AS WELL AS AMBER U.V. PROTECTIVE COATING & AUTOCLAVABLE PLASTIC SAFETY COATING.***

## **PYREX® Brand Proculture® Cell Culture Spinner Flasks:**

- Baffled to enhance agitation of contents of the flask.
- Dynamic Vortex impellers ensure high performance stirring.
- Wide, angled sidearms permit easier entry of 25mL and 50mL pipets for feeding and harvesting.
- Available in 70mm, 100mm, 120mm and 140mm center necks with 32 mm and 45mm sidearms.
- Personalized Flasks -

**We can add Company or School name, lab name, project identifiers, color coding etc.**

**By using the Proculture® Spinner Flask System you are investing in a system that will provide increased efficiency, ease of use, better performance with many cell lines and state of the art cell culture capabilities.**

**Our comprehensive and easy to follow 45 page instruction manual guides you through setup, use and safety procedures. Go to [www.corning.com/lifesciences](http://www.corning.com/lifesciences).**

**Wilmad-LabGlass is the sole supplier of genuine replacement parts for your Corning® Spinner Products.**

## Magnetic Drive Assemblies, Accessories and Components

Impeller Assemblies, Magnetic Drive	CATALOG NUMBER	LIST PRICE
Impeller Assembly, 125mL, Polysulfone Shaft	BP-3401-100	\$67.43
Impeller Assembly, 250mL, Polysulfone Shaft	BP-3401-102	74.36
Impeller Assembly, 500mL, Polysulfone Shaft	BP-3401-104	76.12
Impeller Assembly, 1L, Polysulfone Shaft	BP-3401-106	102.08
Impeller Assembly, 3L, Glass Shaft	BP-3401-108	143.26
Impeller Assembly, 6L, Glass Shaft	BP-3401-110	148.50
Impeller Assembly, 8L, Glass Shaft	BP-3401-112	154.00
Impeller Assembly, 15L, Glass Shaft	BP-3401-114	175.76
Impeller Assembly, 36L, Glass Shaft	BP-3401-116	187.10

Center Neck Cap Assemblies	CATALOG NUMBER	LIST PRICE
Cap Assembly, 70mm, Magnetic, 1/4" Shaft	BP-3402-100	\$66.09
Cap Assembly, 100mm, Magnetic, 1/4" Shaft	BP-3403-100	78.23
Cap Assembly, 100mm, Magnetic, 5/16" Shaft	BP-3404-100	90.16
Cap Assembly, 120mm, Magnetic, 5/16" Shaft	BP-3502-099	102.90
Cap Assembly, 140mm, Magnetic, 5/16" Shaft	BP-3602-100	278.04

## Direct Drive Assemblies, Accessories and Components

Direct Drive Motor	CATALOG NUMBER	LIST PRICE
Motor, 110VAC 60Hz, Complete w/ Connecting Coupling & Drive Yoke	BP-3900-9100	\$1,199.95
Motor, 230VAC 50Hz, Complete w/ Connecting Coupling & Drive Yoke	BP-3900-9102	1,403.12

Cap Assemblies, Direct Drive	CATALOG NUMBER	LIST PRICE
Direct Overhead Drive, 100mm Center Cap, Complete Assembly	BP-3900-9104	\$1,086.87
Direct Overhead Drive, 120mm Center Cap, Complete Assembly	BP-3900-9106	1,190.50

Impeller Assemblies, Direct Drive, Stainless Steel Shaft, w/Coupler	CATALOG NUMBER	LIST PRICE
Impeller Assembly, Direct Drive, 8L	BP-3701-112	\$169.68
Impeller Assembly, Direct Drive, 15L	BP-3701-114	171.59
Impeller Assembly, Direct Drive, 36L	BP-3701-116	195.60

## Proculture® Sidearm Fittings, Angled Sidearm Flasks, Stainless Steel Tubes




500mL Angled Sidearm Flasks	CATALOG NUMBER	LIST PRICE
Gas and Media Handling, 500mL, 1/8" and 1/4" Angled SS Tubes	BP-3830-096	\$220.29

1L Angled Sidearm Flasks	CATALOG NUMBER	LIST PRICE
Gas and Media Handling, 1L, Dual 1/8" Angled Tubes	BP-3824-098	\$199.09
Gas and Media Handling, 1L, Dual 1/4" Angled Tubes	BP-3828-098	199.09
Gas and Media Handling, 1L, Dual 1/8" (1) Long Angled, (1) Straight Tube	BP-3823-098	199.09
Gas and Media Handling, 1L, Dual 1/4" (1) Long Angled, (1) Straight Tube	BP-3827-098	199.09
Gas and Media Handling, 1L, Triple 1/8" (2) Long Angled, (1) Straight Tube	BP-3833-098	215.52
Gas and Media Handling, 1L, Combination (1) 1/8" & (1) 1/4" Angled Tube	BP-3820-098	223.05

3L Angled Sidearm Flasks	CATALOG NUMBER	LIST PRICE
Gas and Media Handling, 3L, Dual 1/8" Angled Tubes	BP-3824-100	\$201.72
Gas and Media Handling, 3L, Dual 1/4" Angled Tubes	BP-3828-100	201.72
Gas and Media Handling, 3L, Dual 1/8" (1) Long Angled, (1) Straight Tube	BP-3823-100	201.72
Gas and Media Handling, 3L, Dual 1/4" (1) Long Angled, (1) Straight Tube	BP-3827-100	201.72
Gas and Media Handling, 3L, Triple 1/8" (2) Long Angled, (1) Straight Tube	BP-3833-100	404.87

6L Angled Sidearm Flasks	CATALOG NUMBER	LIST PRICE
Gas and Media Handling, 6L, Dual 1/8" Angled Tubes	BP-3824-102	\$207.00
Gas and Media Handling, 6L, Dual 1/4" Angled Tubes	BP-3828-102	207.00
Gas and Media Handling, 6L, Dual 1/8" (1) Long Angled, (1) Straight Tube	BP-3823-102	207.00
Gas and Media Handling, 6L, Dual 1/4" (1) Long Angled, (1) Straight Tube	BP-3827-102	207.00
Gas and Media Handling, 6L, Triple 1/8" (2) Long Angled, (1) Straight Tube	BP-3833-102	222.12

8L & 15L Angled Sidearm Flasks	CATALOG NUMBER	LIST PRICE
Gas and Media Handling, 8L, Combination (1) 1/8" & (1) 1/4" Angled Tube	BP-3820-104	\$223.04
Gas and Media Handling, 8L, Dual 1/8" Angled Tubes	BP-3824-104	212.91
Gas and Media Handling, 8L, Dual 1/4" Angled Tubes	BP-3828-104	212.91
Gas and Media Handling, 8L, Dual 1/8" (1) Long Angled, (1) Straight Tube	BP-3823-104	212.91
Gas and Media Handling, 8L, Dual 1/4" (1) Long Angled, (1) Straight Tube	BP-3827-104	212.91
Gas and Media Handling, 8L, Triple 1/8" (2) Long Angled, (1) Straight Tube	BP-3833-104	227.10
Gas and Media Handling, 15L, Dual 1/4" Angled Tubes	BP-3828-106	240.90
Gas and Media Handling, 15L, Dual 1/4" (1) Long Angled, (1) Straight Tube	BP-3827-106	240.90




## Proculture® Sidearm Fittings, Vertical Sidearm Flasks, Stainless Steel Tubes

Gas Handling Fittings	CATALOG NUMBER	LIST PRICE
Gas Delivery Fitting, 1/8" Inlet . . . . .	BP-3801-100 . . . . .	\$163.45
Gas Delivery Fitting, 1/4" Inlet . . . . .	BP-3801-102 . . . . .	172.13
Gas Delivery and Venting Fitting, 1/8" & 1/4" Inlets . . . . .	BP-3801-104 . . . . .	274.74
Gas Venting Cap with 50mm, 0.2µm PTFE Filter . . . . .	BP-3806-102 . . . . .	101.93
Gas Venting Cap with 50mm, 0.2µm PTFE Filter with Sanitary Fitting . . . . .	BP-4519-177 . . . . .	135.86

Insertion Probe Fittings	CATALOG NUMBER	LIST PRICE
12mm Oxygen Probe . . . . .	BP-3811-100 . . . . .	\$109.22
12mm pH Probe . . . . .	BP-3816-100 . . . . .	109.22
12mm Temperature Probe . . . . .	BP-3841-100 . . . . .	109.22
18mm pH, Oxygen or Temperature Probe . . . . .	BP-4519-172 . . . . .	135.79

Media Handling Fittings	CATALOG NUMBER	LIST PRICE
Dual 1/8" Inlets for 8L or 15L Flasks . . . . .	BP-3821-100 . . . . .	\$258.18
Dual 1/8" Inlets for 36L Flasks . . . . .	BP-3821-102 . . . . .	258.18
Dual 1/4" Inlets for 8L or 15L Flasks . . . . .	BP-3826-100 . . . . .	271.44
Dual 1/4" Inlets for 36L Flasks . . . . .	BP-3826-102 . . . . .	271.44
Combination 1/8" and 1/4" for 8L or 15L Flasks . . . . .	BP-3831-100 . . . . .	263.17
Combination 1/8" and 1/4" for 36L Flasks . . . . .	BP-3831-102 . . . . .	263.17
Single 1/2" Inlet for 8L or 15L Flasks . . . . .	BP-3836-100 . . . . .	215.16
Single 1/2" Inlet for 36L Flasks . . . . .	BP-3836-102 . . . . .	231.70
Dual 1/4" Inlets for 8L or 15L Flasks w/ High Temperature EPDM O-ring . . . . .	BP-3826-101 . . . . .	248.52



BP-Vertical Sidearm Fitting

### Replacement Parts for Proculture® Spinner Flasks

Replacement Flasks Flasks Only, 70mm Center Neck, Angled Sidearm	CATALOG NUMBER	LIST PRICE
Flask Only, 125mL Spinner, 70mm, 32mm Sidearms . . . . .	BP-3405-100 . . . . .	\$130.15
Flask Only, 250mL Spinner, 70mm, 32mm Sidearms . . . . .	BP-3405-102 . . . . .	165.33

Flasks Only, 100mm Center Neck, Angled Sidearms	CATALOG NUMBER	LIST PRICE
Flask Only, 500mL Spinner, 100mm, 45mm Sidearms . . . . .	BP-3405-104 . . . . .	\$235.00
Flask Only, 1L Spinner, 100mm, 45mm Sidearms . . . . .	BP-3405-106 . . . . .	290.18
Flask Only, 3L Spinner, 100mm, 45mm Sidearms . . . . .	BP-3405-108 . . . . .	456.43
Flask Only, 6L Spinner, 100mm, 45mm Sidearms . . . . .	BP-3405-110 . . . . .	487.08
Flask Only, 8L Spinner, 100mm, 45mm Sidearms . . . . .	BP-3405-112 . . . . .	551.65
Flask Only, 15L Spinner, 100mm, 45mm Sidearms . . . . .	BP-3405-114 . . . . .	792.66
Flask Only, 36L Spinner, 100mm, 45mm Sidearms . . . . .	BP-3405-116 . . . . .	1,280.62

Flasks Only, 100mm Center Neck, Vertical Sidearms	CATALOG NUMBER	LIST PRICE
Flask Only, 8L Spinner, 100mm, (4) 45mm Sidearms . . . . .	BP-3700-112 . . . . .	\$843.80
Flask Only, 15L Spinner, 100mm, (4) 45mm Sidearms . . . . .	BP-3700-114 . . . . .	1,029.44
Flask Only, 36L Spinner, 100mm, (4) 45mm Sidearms . . . . .	BP-3700-116 . . . . .	1,676.16

Flasks Only, 120mm Center Neck, Angled Sidearms	CATALOG NUMBER	LIST PRICE
Flask Only, 3L Spinner, 120mm, 45mm Sidearms . . . . .	BP-3505-108 . . . . .	\$480.30
Flask Only, 6L Spinner, 120mm, 45mm Sidearms . . . . .	BP-3505-110 . . . . .	507.54
Flask Only, 8L Spinner, 120mm, 45mm Sidearms . . . . .	BP-3505-112 . . . . .	583.88
Flask Only, 15L Spinner, 120mm, 45mm Sidearms . . . . .	BP-3505-114 . . . . .	799.48
Flask Only, 36L Spinner, 120mm, 45mm Sidearms . . . . .	BP-3505-116 . . . . .	1,290.08

Flasks Only, 120mm Center Neck, Vertical Sidearms	CATALOG NUMBER	LIST PRICE
Flask Only, 8L Spinner, 120mm, (4) 45mm Sidearms . . . . .	BP-3700-212 . . . . .	\$865.13
Flask Only, 15L Spinner, 120mm, (4) 45mm Sidearms . . . . .	BP-3700-214 . . . . .	1,063.19
Flask Only, 15L Spinner, 120mm, (6) 45mm Sidearms . . . . .	BP-3700-215 . . . . .	1,255.00
Flask Only, 36L Spinner, 120mm, (6) 45mm Sidearms . . . . .	BP-3700-216 . . . . .	1,712.21

Flasks Only, 140mm Center Neck, Angled Sidearms	CATALOG NUMBER	LIST PRICE
Flask Only, 3L Spinner, 140mm, 45mm Sidearms . . . . .	BP-3605-108 . . . . .	\$469.59
Flask Only, 6L Spinner, 140mm, 45mm Sidearms . . . . .	BP-3605-110 . . . . .	559.57
Flask Only, 8L Spinner, 140mm, 45mm Sidearms . . . . .	BP-3605-112 . . . . .	571.12
Flask Only, 15L Spinner, 140mm, 45mm Sidearms . . . . .	BP-3605-114 . . . . .	850.85
Flask Only, 15L Spinner, 140mm, (4) 45mm Sidearms . . . . .	BP-3700-314 . . . . .	1,096.94
Flask Only, 36L Spinner, 140mm, 45mm Sidearms . . . . .	BP-3605-116 . . . . .	1,348.71
Flask Only, 36L Spinner, 140mm, (6) 45mm Sidearms . . . . .	BP-3700-316 . . . . .	1,653.87

## Replacement Magnetic Impeller Parts

<b>Paddle / Magnet Holder, Lower Rotating, Magnetic Drive</b>	<b>CATALOG NUMBER</b>	<b>LIST PRICE</b>
Bearing / Holder, 125mL.....	BP-3411-010.....	\$20.35
Bearing / Holder, 250mL.....	BP-3411-011.....	20.24
Bearing / Holder, 500mL.....	BP-3411-012.....	28.49
Bearing / Holder, 1L.....	BP-3411-013.....	30.47
Bearing / Holder, 3L.....	BP-3411-014.....	54.45
Bearing / Holder, 6L.....	BP-3411-015.....	68.10
Bearing / Holder, 8L.....	BP-3411-016.....	86.68
Bearing / Holder, 15L.....	BP-3411-017.....	101.05
Bearing / Holder, 36L.....	BP-3411-018.....	126.60



BP-Bearing Holder

<b>Paddle Blade, Magnetic Drive</b>	<b>CATALOG NUMBER</b>	<b>LIST PRICE</b>
Blade, Paddle, 125mL.....	BP-3411-020.....	\$22.00
Blade, Paddle, 250mL.....	BP-3411-021.....	22.88
Blade, Paddle, 500mL.....	BP-3411-022.....	24.20
Blade, Paddle, 1L.....	BP-3411-023.....	36.30
Blade, Paddle, 3L.....	BP-3411-024.....	38.94
Blade, Paddle, 6L.....	BP-3411-025.....	39.82
Blade, Paddle, 8L.....	BP-3411-026.....	47.80
Blade, Paddle, 15L.....	BP-3411-027.....	55.50
Blade, Paddle, 36L.....	BP-3411-028.....	70.10



BP-Paddle Blade

<b>Thrust Washer, Shaft Retaining, Magnetic Drive</b>	<b>CATALOG NUMBER</b>	<b>LIST PRICE</b>
Thrust Washer, 1/4" Diameter, 125mL - 1L Impellers.....	BP-3411-030.....	\$2.09
Thrust Washer, 5/16" Diameter, 3L - 36L Impellers.....	BP-3411-031.....	3.52



BP-Thrust Washer

<b>Paddle Retaining Sleeves</b>	<b>CATALOG NUMBER</b>	<b>LIST PRICE</b>
Retaining Sleeve, 125mL, 250mL, & 500mL.....	BP-3411-050.....	\$11.88
Retaining Sleeve, 1L.....	BP-3411-053.....	13.53
Retaining Sleeve, 3L, 6L, & 8L.....	BP-3411-055.....	18.81
Retaining Sleeve, 15L & 36L.....	BP-3411-057.....	25.41



BP-Retaining Sleeve

<b>Magnets, Impeller Drive, PTFE Encapsulated</b>	<b>CATALOG NUMBER</b>	<b>LIST PRICE</b>
Magnet, 125mL.....	BP-3411-068.....	\$6.38
Magnet, 250mL & 500mL.....	BP-3411-060.....	12.65
Magnet, 1L.....	BP-3411-061.....	13.64
Magnet, 3L - 15L.....	BP-3411-062.....	34.21
Magnet, 36L.....	BP-3411-063.....	60.17



BP-Magnet

<b>Plug, Impeller Sealing, Silicone</b>	<b>CATALOG NUMBER</b>	<b>LIST PRICE</b>
Plug, 125mL - 1L.....	BP-3411-090.....	\$3.27
Plug, 3L - 36L.....	BP-3411-095.....	4.73



BP-Plug

<b>Shaft, Impeller Suspending, Polysulfone</b>	<b>CATALOG NUMBER</b>	<b>LIST PRICE</b>
Shaft, 125mL & 250mL, w/ Thrust Washers.....	BP-3411-170.....	\$9.20
Shaft, 500mL, w/ Thrust Washers.....	BP-3411-172.....	12.20
Shaft, 1L, w/ Thrust Washers.....	BP-3411-173.....	13.00



BP-Polysulfone Shaft

<b>Shaft, Impeller Suspending, Glass</b>	<b>CATALOG NUMBER</b>	<b>LIST PRICE</b>
Shaft, 3L.....	BP-3411-174.....	\$13.29
Shaft, 6L.....	BP-3411-175.....	20.95
Shaft, 8L & 36L.....	BP-3411-176.....	27.80
Shaft, 15L.....	BP-3411-177.....	31.75



BP-Glass Shaft

<b>Shaft, Impeller Suspending, Stainless Steel</b>	<b>CATALOG NUMBER</b>	<b>LIST PRICE</b>
Shaft, 1L.....	BP-3611-083.....	\$49.17
Shaft, 3L.....	BP-3611-084.....	59.90
Shaft, 6L.....	BP-3611-085.....	64.70
Shaft, 8L.....	BP-3611-086.....	72.00
Shaft, 15L.....	BP-3611-087.....	88.80
Shaft, 36L.....	BP-3611-088.....	101.40



BP-Stainless Steel Shaft

## Replacement Center Neck Parts

70mm Center Neck Cap	CATALOG NUMBER	LIST PRICE
70mm Cap, Screw Type, w/ Integrated Compression Fitting (ICF) . . . . .	BP-3412-011 . . . . .	\$39.95
Liner, 70mm w/ Center Hole, Red Silicone . . . . .	BP-3412-015 . . . . .	14.50



BP-CAP-ICF

100mm Center Neck Cap	CATALOG NUMBER	LIST PRICE
100mm Cap, Screw Type, w/ Integrated Compression Fitting (ICF) . . . . .	BP-3413-011 . . . . .	\$56.10
Liner, 100mm w/ Center Hole, Red Silicone . . . . .	BP-3413-015 . . . . .	18.50

120mm Center Neck Cap	CATALOG NUMBER	LIST PRICE
120mm Cap, Screw Type, w/ Integrated Compression Fitting (ICF) . . . . .	BP-3516-011 . . . . .	\$70.00
Liner, 120mm w/ Center Hole, Red Silicone . . . . .	BP-3513-015 . . . . .	20.50

140mm Center Neck Cap	CATALOG NUMBER	LIST PRICE
140mm Cap, Lower Plate . . . . .	BP-3612-010 . . . . .	\$150.05
140mm Cap, Upper Plate . . . . .	BP-3612-020 . . . . .	95.57
140mm Cap, Seal Actuating Ring . . . . .	BP-3612-030 . . . . .	39.82
140mm Cap, Washer, Actuating Ring Retaining . . . . .	BP-3612-050 . . . . .	5.06
140mm Cap, O-ring, Actuating Ring Retaining . . . . .	BP-3612-060 . . . . .	3.41
140mm Cap, Gasket, Compression Sealing . . . . .	BP-3612-070 . . . . .	12.98



BP-LINER

Nut, Impeller Retaining, Ultem™	CATALOG NUMBER	LIST PRICE
Nut, 1/4", 125mL - 1L . . . . .	BP-3413-9039 . . . . .	\$18.59
Nut, 5/16", 3L - 36L . . . . .	BP-3413-9040 . . . . .	18.59



BP-NUT

Autoclavable Glass Reinforced PBT Storage Caps w/Silicone Liner	CATALOG NUMBER	LIST PRICE
70mm Cap Assembly . . . . .	BP-3412-017 . . . . .	\$23.75
100mm Cap Assembly . . . . .	BP-3413-018 . . . . .	26.25
120mm Cap Assembly . . . . .	BP-3513-018 . . . . .	51.60



BP-Storage Cap

### Replacement Parts for Proculture® Sidearm Fittings

PET Inserts	CATALOG NUMBER	LIST PRICE
45mm PET Insert, Single 1/8" Inlet . . . . .	BP-3802-001 . . . . .	\$32.90
45mm PET Insert Single 1/4" Inlet . . . . .	BP-3803-001 . . . . .	33.12
45mm PET Insert for 50mm Filter . . . . .	BP-3807-001 . . . . .	53.51
45mm PET Insert for 12mm Insertion Probes . . . . .	BP-3812-001 . . . . .	53.51
45mm PET Insert, Dual 1/8" Inlet . . . . .	BP-3822-001 . . . . .	53.51
45mm PET Insert, Dual 1/4" Inlet . . . . .	BP-3827-001 . . . . .	53.51
45mm PET Insert, Combination 1/8" and 1/4" Inlet . . . . .	BP-3832-001 . . . . .	58.51
45mm PET Insert, Triple 1/8" Inlet . . . . .	BP-3832-002 . . . . .	63.10



BP-Insert

Tubes, Type 316 Stainless Steel For Vertical Sidearm Flasks	CATALOG NUMBER	LIST PRICE
1/8" O.D. X 18" Long For 8L and 15L Flasks . . . . .	BP-3802-007 . . . . .	\$44.15
1/4" O.D. X 18" Long w/ 1/4"-28 Thread for 8L and 15L Flasks . . . . .	BP-3803-007 . . . . .	77.70
1/8" O.D. X 21" Long For 36L Flasks . . . . .	BP-3822-007 . . . . .	56.75
1/4" O.D. X 21" Long For 36L Flasks . . . . .	BP-3827-007 . . . . .	69.80

*For Replacement Angled Sidearm Flasks Tubes Contact Customer Service.*


### Misc Replacement Parts & Fittings

Compression Nuts and Port Rings	CATALOG NUMBER	LIST PRICE
Nut, 1/8", Noryl, w/ Integrated Ferrule . . . . .	BP-3802-004 . . . . .	\$12.42
Nut, 1/4", Noryl, w/ Integrated Ferrule . . . . .	BP-3803-004 . . . . .	12.42
Nut, 12mm, Noryl . . . . .	BP-3812-004 . . . . .	16.21
Plug, 1/8", Acetal . . . . .	BP-3847-008 . . . . .	19.46
Plug, 1/4", Acetal . . . . .	BP-3852-008 . . . . .	20.75
Plug, 1/2" & 12mm, Acetal . . . . .	BP-3857-008 . . . . .	23.05



BP-Nut Noryl

Caps, O-rings, Filters, and Misc. Parts	CATALOG NUMBER	LIST PRICE
Cap, 45mm, Sidearm Securing, Orange . . . . .	BP-3802-003 . . . . .	\$10.65
Cap, 45mm, Solid Top, Orange . . . . .	C-1395-45LTC . . . . .	2.24
Cap, 32mm, Solid Top, Orange . . . . .	C-1395-32LTC . . . . .	2.40
Drip Ring, 45mm, Clear . . . . .	C-1395-45LTR . . . . .	1.02
Sparging Stone & Nut, Stainless Steel, 1/8" Tubing or 1/4"-28 Thread . . . . .	BP-3802-006 . . . . .	49.50
Filter, 50mm, 2µm, PTFE . . . . .	BP-3807-004 . . . . .	21.08
O-ring, Viton, Insert Sealing, 026 . . . . .	BP-3802-002 . . . . .	2.00
O-ring, Viton, 12mm Compression Fitting . . . . .	BP-3812-005 . . . . .	1.50



BP-Cap 45mm